

YZ Pump Head

YZ15-13X/YZ25-13X



YZ15-13X/YZ25-13X vs. YZ1515x/YZ2515x



Main advantages of YZ15-13X/YZ25-13X

- 1. The extended roller design supports larger tubing sizes for improved tubing compatibility and reliable operation.
- 2. Automatic tuing clamp lifting for easier tubing loading.
- 3. The tubing clamp conforms gently to the tubing, reducing stress and minimizing the risk of rupture.

Model	Prpduct Code	Roller Number	Roller Material	Housing Marterial	Speed (rpm)	Tubing	Max. Flow Rate (mL/min)	Weight (kg)
YZ15-13X	05.01.54A	3	Stainless steel	PESU	13# 14# 19# 16# 25# 17# 18#	≤650rpm	2200mL/min @600rpm 2350mL/min@650rpm	0.46
YZ25-13X	05.01.57A	3	Stainless steel	PESU	15# 24# 35# 36#	≤650rpm	3000mL/min @600rpm 3300mL/min@650rpm	0.46



- Reliable classic pump head design compatible with multiple tubing sizes.
- Wide flow range with simple and efficient operation.

